



Mob3A/B Monoclonal Antibody

Catalog No	BYmab-04345
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	MMOBKL2A/MMOBKL2B
Protein Name	Mps one binder kinase activator-like 2A/B
Immunogen	Synthesized peptide derived from Mob3A/B . at AA range: 1-80
Specificity	Mob3A/B Monoclonal Antibody detects endogenous levels of Mob3A/B protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MOB3A; MOBKL2A; MOB kinase activator 3A; MOB-LAK; Mob1 homolog 2A; Mps one binder kinase activator-like 2A; MOB3B; MOBKL2B; MOB kinase activator 3B; Mob1 homolog 2b; Mps one binder kinase activator-like 2B
Synonyms Observed Band	Mps one binder kinase activator-like 2A; MOB3B; MOBKL2B; MOB kinase
	Mps one binder kinase activator-like 2A; MOB3B; MOBKL2B; MOB kinase activator 3B; Mob1 homolog 2b; Mps one binder kinase activator-like 2B
Observed Band	Mps one binder kinase activator-like 2A; MOB3B; MOBKL2B; MOB kinase activator 3B; Mob1 homolog 2b; Mps one binder kinase activator-like 2B 25kD
Observed Band Cell Pathway	Mps one binder kinase activator-like 2A; MOB3B; MOBKL2B; MOB kinase activator 3B; Mob1 homolog 2b; Mps one binder kinase activator-like 2B 25kD intracellular,
Observed Band Cell Pathway Tissue Specificity	Mps one binder kinase activator-like 2A; MOB3B; MOBKL2B; MOB kinase activator 3B; Mob1 homolog 2b; Mps one binder kinase activator-like 2B 25kD intracellular, Cervix,Skin, function:May regulate the activity of kinases.,sequence caution:Chimera. There is an inclusion of DNA from another gene in position 141, similarity:Belongs to the

Nanjing BYabscience technology Co.,Ltd



国内优质抗体供应商 精准的 WB 检测服务 24H 在线服务,欢迎咨询



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

