



KOR-1 (phospho Ser369) Monoclonal Antibody

Catalog No	BYmab-12638
lsotype	lgG
Reactivity	Mouse;Rat
Applications	WB
Gene Name	OPRK1
Protein Name	Kappa-type opioid receptor
Immunogen	The antiserum was produced against synthesized peptide derived from mouse KOR-1 around the phosphorylation site of Ser369. AA range:331-380
Specificity	Phospho-KOR-1 (S369) Monoclonal Antibody detects endogenous levels of KOR-1 protein only when phosphorylated at S369.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	OPRK1; OPRK; Kappa-type opioid receptor; K-OR-1; KOR-1
Observed Band	42kD
Cell Pathway	
Tissue Specificity	
Function	
Background	Endogenous opioid peptides and opiates, like morphine, transmit their pharmacological effects through membrane bound opioid receptors. Pharmacological studies and molecular cloning have led to the identification of three different types of opioid receptor, mu-type, delta-type and kappa-type, also designated MOR-1, DOR-1 and KOR-1, respectively. MOR-1 is a receptor for beta-endorphin, DOR-1 is a receptor for enkephalins, and KOR-1 is a receptor for dynorphins.The three opioid receptor types are highly homologous and belong to the superfamily of G-protein-coupled receptors. Opioid receptors have been

Nanjing BYabscience technology Co.,Ltd




国内优质抗体供应商

国内1/L质加神铁应的 — 24H 在线服务, 欢迎咨询 — 24H 在线服务, 欢迎咨询 — 24H 在线服务, 欢迎咨询

shown to modulate a range of brain functions, including instinctive behavior and
emotions. This regulation is thought to involve the inhibition of neurotransmitter
release by reducing calcium ion currents and increasing potassium ion
conductance.matters needing
attentionAvoid repeated freezing and thawing!Usage suggestionsThis product can be used in immunological reaction related experiments. For
more information, please consult technical personnel.

Products Images



Nanjing BYabscience technology Co.,Ltd